inno ge

oxygen. anytime. anywhere.™

**MEDITEK**  
Medisinsk utstyr

# Oksygenkonsentratoren som definerer en ny utstyrskategori

---

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# Titration Summary 2011 Inogen One G2

## Titration Summary 2011

- Represents 5396 titrations (5396 different patients)
- All titrations completed with either licensed Respiratory Therapist or Nurse
- Target patients 3 LPM or less prescription
- Patient must achieve 90% or 92% oxygen saturation (specified by doctor) at rest and in activity, OR equivalent saturation compared to old system to be considered “pass/successful”
- Whatever flow setting is found to achieve saturation at activity levels, is designate as the nocturnal setting. Past trials indicated that the activity setting is the appropriate setting to be used at night<sup>1</sup>

<sup>1</sup> (“Determination of an appropriate nocturnal setting for a portable oxygen concentrator with a pulse-dosed delivery” revA Stegmaier Inogen Nocturnal Study – 2006 Final).

# Overall figures 2011 – Inogen One G2

Number of Inogen One G2 titrations completed	Titration success rate for Inogen One G2 %
Jan: 333	99,4 %
Feb: 391	99,7 %
Mar: 508	99,8 %
April: 363	99,7 %
May: 335	99,7 %
June: 358	99,2 %
July: 465	98,9 %
Aug: 552	99,6 %
Sep: 470	99,8 %
Oct: 595	98,2 %
Nov: 496	99,2 %
Dec: 530	98,5 %
<b>Total titrations 2011: 5396</b>	<b>Total success rate % 2011: 99,30</b>

## Product that has defined a new category

---

1. Stegmaier JP, Chatburn RL, Lewarski JS. "Determination of an Appropriate Nocturnal Setting for a Portable Oxygen Concentrator with Pulsed-Dosed Delivery." Abstract. *Respir Care* November 2006;51(11): 1305.

Summary: The purpose of this study was to determine if a single titration of oxygen using a POC during ambulation/exercise would provide an appropriate setting for nocturnal use. The results suggested that an oxygen setting selection based on daytime ADL/ambulation appears to produce effective nocturnal oxygen therapy as evidenced by a mean sleeping SpO<sub>2</sub> of 92% and no clinically significant desaturation.

2. Chatburn, R, Lewarski J, McCoy R. "Nocturnal oxygenation using a pulsed dose oxygen conserving device compared to continuous flow oxygen." *Respir Care* March 2006;51(3): 252-256

Summary: The study compared nocturnal oxygenation with continuous flow versus the Inogen One among a group of established LTOT users with chronic lung disease. The results demonstrate that when appropriately titrated, the Inogen One is essentially clinically equivalent to continuous flow oxygen. The study also suggests that daytime pulse dose titrations may be effective in determining nocturnal oxygenation.

## Product that has *defined* a new category

---

3. Case, R, Hausmann R. "Use of a portable oxygen concentrator with a fixed minute volume oxygen conserving device to deliver oxygen to exercising pulmonary rehabilitation patients." Abstract.

Respir Care November 2005;50(11):1510.

Summary: The study concluded that the Inogen One was as clinically effective as continuous flow oxygen at maintaining target SpO<sub>2</sub> levels in high flow (4-5 lpm) oxygen users during intense exercise.

4. McCoy, R, Lewarski, J. "A test for clinical equivalency: A portable concentrator with integrated oxygen-conserving compared to continuous flow oxygen during nocturnal use." Abstract.

Respir Care November 2005; 51 (11):

Summary: This study demonstrated the Inogen One™ POC was able to deliver adequate nocturnal oxygen therapy as evidenced by continuous SpO<sub>2</sub> monitoring in 9 of 10 (90%) of patients studied. The resting daytime oxygen titration and the resultant SpO<sub>2</sub> appears to be reasonably effective method for determining an appropriate nocturnal oxygen setting.

# Awards 2011

**“Industrial designers Society of America”**

**Category: “Solution to a developed world social problem”**

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